

Anti-CD49b/pan-NK cells antibody [DX5] {EV450} (STJA0016406)

STJA0016406

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin Alpha-2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DX5
Concentration	5 Mu L/Test
Conjugation	EV450
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	16398
Gene Symbol	Itga2
Uniprot ID	ITA2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	