

Anti-CD49b/pan-NK cells antibody [DX5] {ER780} (STJA0016407)

STJA0016407

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Integrin Alpha-2 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | DX5 |
| Concentration | 5 µg L/Test |
| Conjugation | ER780 |
| Purification | Affinity purified |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgMk |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 16398 |
| Gene Symbol | Itga2 |
| Uniprot ID | ITA2_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |