

Anti-Ly-6G/Ly-6C (Gr-1) antibody [RB6-8C5] {EF700} (STJA0016410)

STJA0016410

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ly-6G/Ly-6C (Gr-1) is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RB6-8C5
Concentration	
Conjugation	EF700
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	100041546
	546644
Gene Symbol	Ly6c2
	Ly6g
Uniprot ID	LY6C2_MOUSE
	LY6G_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	