

Anti-Ly6C antibody [Monts 1] {EV610} (STJA0016412)

STJA0016412

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Lymphocyte antigen 6c2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Monts 1
Concentration	
Conjugation	EV610
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	100041546
Gene Symbol	Ly6c2
Uniprot ID	LY6C2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	