

Anti-HLA-A, B, C antibody [W6/32] {EV450} (STJA0016418)

STJA0016418

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-A, B, C is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	W6/32
Concentration	5 µg L/Test
Conjugation	EV450
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	3106 3107 3105 HLA-B HLA-C HLAB_HUMAN HLAC_HUMAN HLAA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	