

Anti-CD274/PD-L1 antibody [10F.9G2] {EF488} (STJA0016420)

STJA0016420

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Programmed Cell Death 1 Ligand 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10F.9G2
Concentration	5 Mu L/Test
Conjugation	EF488
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	60533
Gene Symbol	Cd274
Uniprot ID	PD1L1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	