

Anti-CD45 antibody [30-F11] {EF700} (STJA0016421)

STJA0016421

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Receptor-Type Tyrosine-Protein Phosphatase C is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	30-F11
Concentration	
Conjugation	EF700
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	19264
Gene Symbol	Ptpnc
Uniprot ID	PTPRC_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	