

## Anti-CD150/SLAM antibody [TC15-12F12.2] {EF647} (STJA0016443)

STJA0016443

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Signaling Lymphocytic Activation Molecule is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	TC15-12F12.2
<b>Concentration</b>	5 Mu L/Test
<b>Conjugation</b>	EF647
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG2aL
<b>Storage Instruction</b>	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">27218</a>
<b>Gene Symbol</b>	<a href="#">Slamf1</a>
<b>Uniprot ID</b>	<a href="#">SLAF1_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	