

Anti-CD23 antibody [B3B4] {EF647} (STJA0016446)

STJA0016446

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Low Affinity Immunoglobulin Epsilon Fc Receptor is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B3B4
Concentration	5 Mu L/Test
Conjugation	EF647
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	14128
Gene Symbol	Fcer2
Uniprot ID	FCER2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	