

Anti-IgD antibody [11-26c.2a] {EV450} (STJA0016459)

STJA0016459

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ig Delta Chain C Region Secreted Form is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11-26c.2a
Concentration	5 Mu L/Test
Conjugation	EV450
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	IGHD_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	