

Anti-Ly6A/E (Sca-1) antibody [D7] {EV450} (STJA0016462)

STJA0016462

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Lymphocyte antigen 6a-2/6e-1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	D7
Concentration	5 Mu L/Test
Conjugation	EV450
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	110454
Gene Symbol	Ly6a
Uniprot ID	LY6A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	