

Anti-H-2 antibody [M1/42] {EV450} (STJA0016472)

STJA0016472

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-H-2 Class II Histocompatibility antigen-A-S Beta Chain is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M1/42
Concentration	5 Mu L/Test
Conjugation	EV450
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	
Gene Symbol	H2-Ab1
Uniprot ID	HB2S_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	