

Anti-CD200/OX2 antibody [OX-90] {FITC} (STJA0016479)

STJA0016479

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ox-2 Membrane Glycoprotein is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-90
Concentration	5 Mu L/Test
Conjugation	FITC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	
Gene Symbol	Cd200
Uniprot ID	OX2G_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	