

Anti-CD16 antibody [3G8] {EV610} (STJA0016481)

STJA0016481

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD16 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3G8
Concentration	
Conjugation	EV610
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	2214 2215
Gene Symbol	FCGR3A FCGR3B
Uniprot ID	FCG3A_HUMAN FCG3B_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	