

Anti-CD8a antibody [HIT8a] {FITC} (STJA0016492)

STJA0016492

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-Cell Surface Glycoprotein Cd8 Alpha Chain is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HIT8a
Concentration	5 Mu L/Test
Conjugation	FITC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	925
Gene Symbol	CD8A
Uniprot ID	CD8A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	