

Anti-TCR Gamma/Delta antibody [GL3] {FITC} (STJA0016496)

STJA0016496

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-TCR Gamma/Delta is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GL3
Concentration	5 µg L/Test
Conjugation	FITC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	10219
	50928
Gene Symbol	KLRG1
	Klrg1
Uniprot ID	KLRG1_HUMAN
	KLRG1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	