

Anti-CD172a/SIRP Alpha antibody [P84] {EF488} (STJA0016498)

STJA0016498

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD172a/SIRP Alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	P84
Concentration	
Conjugation	EF488
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	12490
Gene Symbol	Cd34
Uniprot ID	CD34_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	