

Anti-CD163 antibody [S15049F] {APC} (STJA0016500)

STJA0016500

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Scavenger Receptor Cysteine-Rich Type 1 Protein M130 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S15049F
Concentration	5 Mu L/Test
Conjugation	APC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	93671
Gene Symbol	Cd163
Uniprot ID	C163A_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	