

Anti-CD203c antibody [NP4D6] {EF488} (STJA0016503)

STJA0016503

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Ectonucleotide Pyrophosphatase/Phosphodiesterase Family Member 3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	NP4D6
Concentration	5 Mu L/Test
Conjugation	EF488
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	5169
Gene Symbol	ENPP3
Uniprot ID	ENPP3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	