

## Anti-CD56/NCAM antibody [B-A19] {PE/Cyanine 7} (STJA0016524)

STJA0016524

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-C-C Chemokine Receptor Type 3 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	B-A19
<b>Concentration</b>	5 Mu L/Test
<b>Conjugation</b>	PE/Cyanine 7
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">1232</a>
<b>Gene Symbol</b>	<a href="#">CCR3</a>
<b>Uniprot ID</b>	<a href="#">CCR3_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	