

Anti-CD45RA antibody [OX-33] {PE} (STJA0016530)

STJA0016530

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD45RA is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | OX-33 |
| Concentration | 5 µg L/Test |
| Conjugation | PE |
| Purification | Affinity purified |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG1k |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze. |

TARGET INFORMATION

| | |
|---------------------------|---------------------------|
| Gene ID | 24699 |
| Gene Symbol | Ptprc |
| Uniprot ID | PTPRC_RAT |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |