

Anti-CD45RA antibody [OX-33] {APC} (STJA0016531)

STJA0016531

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD45RA is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-33
Concentration	5 µg L/Test
Conjugation	APC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	24699
Gene Symbol	Ptprc
Uniprot ID	PTPRC_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	