

Anti-CD161 antibody [3.2.3] {FITC} (STJA0016533)

STJA0016533

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD161 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3.2.3
Concentration	5 Mu L/Test
Conjugation	FITC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	NA 25192 678513 362443
Gene Symbol	Klrb1b Klrb1b KRBBB_RAT KRBBB_RAT KRBBB_RAT <
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	