

**Anti-CD5 antibody [HISM2] {EV450} (STJA0016551)**

STJA0016551

**GENERAL INFORMATION**

|                          |                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                          |
| <b>Short Description</b> | Mouse monoclonal antibody anti-T-Cell Surface Glycoprotein Cd5 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                                          |
| <b>Host/Source</b>       | Mouse                                                                                                                       |
| <b>Reactivity</b>        | Human                                                                                                                       |

**PRODUCT PROPERTIES**

|                            |                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                              |
| <b>Clone ID</b>            | HISM2                                                                                                                   |
| <b>Concentration</b>       |                                                                                                                         |
| <b>Conjugation</b>         | EV450                                                                                                                   |
| <b>Purification</b>        | Affinity purified                                                                                                       |
| <b>Dilution Range</b>      |                                                                                                                         |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                             |
| <b>Isotype</b>             | IgG1k                                                                                                                   |
| <b>Storage Instruction</b> | This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze. |

**TARGET INFORMATION**

|                         |           |
|-------------------------|-----------|
| <b>Gene ID</b>          | 921       |
| <b>Gene Symbol</b>      | CD5       |
| <b>Uniprot ID</b>       | CD5_HUMAN |
| <b>Immunogen</b>        |           |
| <b>Immunogen Region</b> |           |
| <b>Specificity</b>      |           |
| <b>Immunogen</b>        |           |
| <b>Sequence</b>         |           |