

Anti-CD123 antibody [HI12H7] {PE} (STJA0016614)

STJA0016614

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Interleukin-3 Receptor Subunit Alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HI12H7
Concentration	5 Mu L/Test
Conjugation	PE
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	3563
Gene Symbol	IL3RA
Uniprot ID	IL3RA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	