

Anti-Foxp3 antibody [FJK-16s] {APC} (STJA0016620)

STJA0016620

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Forkhead Box Protein P3 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | FJK-16s |
| Concentration | |
| Conjugation | APC |
| Purification | Affinity purified |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze. |

TARGET INFORMATION

| | |
|---------------------------|-----------------------------|
| Gene ID | 20371 |
| Gene Symbol | Foxp3 |
| Uniprot ID | FOXP3_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |