

Anti-Foxp3 antibody [FJK-16s] {PE/EF594} (STJA0016623)

STJA0016623

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Forkhead Box Protein P3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FJK-16s
Concentration	
Conjugation	PE/EF594
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	20371
Gene Symbol	Foxp3
Uniprot ID	FOXP3_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	