

Anti-CD49d antibody [HP1/2] {PerCP/Cyanine 5.5} (STJA0016647)

STJA0016647

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Integrin Alpha-4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HP1/2
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine 5.5
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	3676
Gene Symbol	ITGA4
Uniprot ID	ITA4_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	