

Anti-CD101 antibody [BB27] {PerCP/Cyanine 5.5} (STJA0016651)

STJA0016651

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Immunoglobulin Superfamily Member 2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BB27
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine 5.5
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Do not freeze.

TARGET INFORMATION

Gene ID	9398
Gene Symbol	CD101
Uniprot ID	IGSF2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	