

**Anti-CD194/CCR4 antibody [L291H4] {PerCP/Cyanine 5.5} (STJA0016654)**

STJA0016654

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-C-C Chemokine Receptor Type 4 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	L291H4
<b>Concentration</b>	5 Mu L/Test
<b>Conjugation</b>	PerCP/Cyanine 5.5
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	This product can be stored at 2-8°C for 12 months. Do not freeze.

**TARGET INFORMATION**

<b>Gene ID</b>	1233
<b>Gene Symbol</b>	CCR4
<b>Uniprot ID</b>	CCR4_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	