

Anti-CD226/DNAM-1 antibody [11A8] {PE} (STJA0016656)
STJA0016656**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cd226 antigen is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11A8
Concentration	5 µg L/Tests
Conjugation	PE
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	10666
Gene Symbol	CD226
Uniprot ID	CD226_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	