

Anti-CD244 antibody [C1.7] {APC} (STJA0016674)

STJA0016674

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Natural Killer Cell Receptor 2b4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	C1.7
Concentration	5 Mu L/Tests
Conjugation	APC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	51744
Gene Symbol	CD244
Uniprot ID	CD244_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	