

Anti-CD95/Fas antibody [HFE7A/APO] {PE} (STJA0016685)

STJA0016685

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse recombinant antibody anti-CD95/Fas is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Recombinant |
| Clone ID | HFE7A/APO |
| Concentration | |
| Conjugation | PE |
| Purification | Affinity purified |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 355 |
| Gene Symbol | FAS |
| Uniprot ID | TNR6_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081