

Anti-CD95/Fas antibody [HFE7A/APO] {EF647} (STJA0016690)

STJA0016690

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse recombinant antibody anti-Tumor Necrosis Factor Receptor Superfamily Member 6 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Recombinant
Clone ID	HFE7A/APO
Concentration	5 Mu L/Tests
Conjugation	EF647
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	355
Gene Symbol	FAS
Uniprot ID	TNR6_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	