

## Anti-CD158b/j antibody [DX27] {PerCP/Cyanine 5.5} (STJA0016694)

STJA0016694

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD158b/j is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	DX27
<b>Concentration</b>	5 Mu L/Tests
<b>Conjugation</b>	PerCP/Cyanine 5.5
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG2ak
<b>Storage Instruction</b>	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3803</a> <a href="#">3804</a>
<b>Gene Symbol</b>	<a href="#">KIR2DL2</a> <a href="#">KIR2DL3</a>
<b>Uniprot ID</b>	<a href="#">K12L2_HUMAN</a> <a href="#">K12L3_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	