

Anti-CD23 antibody [EBVCS2] {EV450} (STJA0016700)

STJA0016700

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Low Affinity Immunoglobulin Epsilon Fc Receptor is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBVCS2
Concentration	5 Mu L/Tests
Conjugation	EV450
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	2208
Gene Symbol	FCER2
Uniprot ID	FCER2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	