

Anti-CD94 antibody [DX22] {PE/Cyanine 7} (STJA0016704)

STJA0016704

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Natural Killer Cells antigen Cd94 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DX22
Concentration	5 Mu L/Tests
Conjugation	PE/Cyanine 7
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	3824
Gene Symbol	KLRD1
Uniprot ID	KLRD1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	