

Anti-CD192/CCR2 antibody [K036C2] {APC} (STJA0016708)
STJA0016708**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-C-C Chemokine Receptor Type 2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	K036C2
Concentration	
Conjugation	APC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	729230
Gene Symbol	CCR2
Uniprot ID	CCR2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	