

Anti-CD158/KIR2DL1 antibody [HP-MA4] {APC} (STJA0016712)

STJA0016712

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD158/KIR2DL1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HP-MA4
Concentration	5 µg L/Tests
Conjugation	APC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	3802 3810 3806 KIR2DL1 KIR2DS5 KI2L1_HUMAN KI2S5_HUMAN KI2S1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	