

**Anti-CD158/KIR2DL1 antibody [HP-MA4] {PerCP/Cyanine 5.5} (STJA0016713)**

STJA0016713

**GENERAL INFORMATION**

|                          |                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                        |
| <b>Short Description</b> | Mouse monoclonal antibody anti-CD158/KIR2DL1 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                        |
| <b>Host/Source</b>       | Mouse                                                                                                     |
| <b>Reactivity</b>        | Human                                                                                                     |

**PRODUCT PROPERTIES**

|                            |                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                              |
| <b>Clone ID</b>            | HP-MA4                                                                                                                  |
| <b>Concentration</b>       | 5 Mu L/Tests                                                                                                            |
| <b>Conjugation</b>         | PerCP/Cyanine 5.5                                                                                                       |
| <b>Purification</b>        | Affinity purified                                                                                                       |
| <b>Dilution Range</b>      |                                                                                                                         |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                             |
| <b>Isotype</b>             | IgG2bk                                                                                                                  |
| <b>Storage Instruction</b> | This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze. |

**TARGET INFORMATION**

|                           |                                                                                                                                                                                                                         |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">3802</a><br><a href="#">3806</a><br><a href="#">3810</a><br><a href="#">KIR2DL1</a><br><a href="#">KIR2DS1</a><br><a href="#">KI2L1_HUMAN</a><br><a href="#">KI2S1_HUMAN</a><br><a href="#">KI2S5_HUMAN</a> |
| <b>Immunogen</b>          |                                                                                                                                                                                                                         |
| <b>Immunogen Region</b>   |                                                                                                                                                                                                                         |
| <b>Specificity</b>        |                                                                                                                                                                                                                         |
| <b>Immunogen Sequence</b> |                                                                                                                                                                                                                         |