

Anti-CD226 antibody [10E5] {APC} (STJA0016716)

STJA0016716

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Cd226 antigen is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10E5
Concentration	5 Mu L/Tests
Conjugation	APC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	225825
Gene Symbol	Cd226
Uniprot ID	CD226_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	