

Anti-CD20 antibody [SA271G2] {PE/Cyanine 5} (STJA0016719)

STJA0016719

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-B-Lymphocyte antigen Cd20 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SA271G2
Concentration	5 µg L/Tests
Conjugation	PE/Cyanine 5
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	12482
Gene Symbol	Ms4a1
Uniprot ID	CD20_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	