

Anti-CD185 antibody [MU5UBEE] (PE) (STJA0016738)

STJA0016738

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-C-X-C Chemokine Receptor Type 5 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MU5UBEE
Concentration	5 Mu L/Tests
Conjugation	PE
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	643
Gene Symbol	CXCR5
Uniprot ID	CXCR5_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	