

**Anti-CD146 antibody [AA98] {FITC} (STJA0016741)**

STJA0016741

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Cell Surface Glycoprotein Muc18 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	AA98
<b>Concentration</b>	5 Mu L/Tests
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG2ak
<b>Storage Instruction</b>	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">4162</a>
<b>Gene Symbol</b>	<a href="#">MCAM</a>
<b>Uniprot ID</b>	<a href="#">MUC18_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	