

Anti-CD146 antibody [AA98] {PE/Cyanine 7} (STJA0016745)

STJA0016745

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cell Surface Glycoprotein Muc18 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AA98
Concentration	5 Mu L/Tests
Conjugation	PE/Cyanine 7
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	4162
Gene Symbol	MCAM
Uniprot ID	MUC18_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	