

Anti-CD74 antibody [LL1] {APC} (STJA0016755)

STJA0016755

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Hla Class II Histocompatibility antigen Gamma Chain is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LL1
Concentration	5 µg L/Tests
Conjugation	APC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	972
Gene Symbol	CD74
Uniprot ID	HG2A_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	