

**Anti-CD8a/CD56/CD3/CD45 Cocktail antibody [CB19;5.1H11;OKT-3;HI30]
{FITC/PE/PE/Cyanine 7/PerCP} (STJA0016759)**

STJA0016759

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	NA monoclonal antibody anti-CD8a/CD56/CD3/CD45 Cocktail is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CB19;5.1H11;OKT-3;HI30
Concentration	5 Mu L/Test
Conjugation	FITC/PE/PE/Cyanine 7/PerCP
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence