

## Human FE (FE) Antibody Pair (STJAP00024)

STJAP00024

### GENERAL INFORMATION

<b>Product Type</b>	Antibody Pairs
<b>Short Description</b>	This FE antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rabbit/Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal/Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated/Biotin
<b>Purification</b>	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
<b>Dilution Range</b>	
<b>Formulation</b>	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">2495</a> <a href="#">2512</a>
<b>Gene Symbol</b>	<a href="#">FTH1</a> <a href="#">FTL</a>
<b>Uniprot ID</b>	<a href="#">FRIH_HUMAN</a> <a href="#">FRIL_HUMAN</a>
<b>Immunogen</b>	Capture immunogen: Native Protein. Detection Immunogen: Native Protein.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	