

Human LpPLA2 (Platelet-Activating Factor Acetylhydrolase) Antibody Pair (STJAP00045)

STJAP00045

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This LpPLA2 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7941
Gene Symbol	PLA2G7
Uniprot ID	PAFA_HUMAN
Immunogen	Capture immunogen: Recombinant Human LpPLA2 protein. Detection Immunogen: Recombinant Human LpPLA2 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	