

Human MMP-1 (Interstitial Collagenase) Antibody Pair (STJAP00047)

STJAP00047

GENERAL INFORMATION

Product Type	Antibody Pairs
Short Description	This MMP-1 antibody pair contains complimentary capture and detection antibodies, which work together under the principals of the Sandwich ELISA.
Applications	ELISA
Host/Source	Goat/Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Dilution Range	
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Isotype	
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4312
Gene Symbol	MMP1
Uniprot ID	MMP1_HUMAN
Immunogen	Capture immunogen: Recombinant Human MMP-1 protein. Detection Immunogen: Recombinant Human MMP-1 protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	